

WHEELCHAIR SCALES



Model 350-10-3

**VACUUM FORMED
SCALE TOP**

RICE LAKE®
WEIGHING SYSTEMS

800-472-6703
www.ricelake.com/health

Rice Lake Weighing Systems offers a comprehensive line of safe, durable wheelchair scales designed to provide ultimate support and accuracy. Motiontrap™ movement compensation technology adjusts for involuntary movement and keyboard tare efficiently reduces errors. Built with a wide range of handrail, seat and capacity options, Rice Lake's line of wheelchair scales simplifies patient weighing.



350-10 Series

Wheelchair Scales

Standard Features

- Capacity: 1000 lb x 0.2 lb (450 kg x 0.1 kg)
- Platform weighing dimensions: 29.5 in x 32 in
- Handle and two wheels for ease of portability
- Motiontrap™ (movement compensation technology)
- Large 1-inch LCD display
- Push-button or manual entry of tare weight
- BMI (body mass index) function
- Hold function
- Unit of measure toggle or locking
- Non-skid platform
- EMR capable, USB
- Powered by 6 AA batteries (included)

Options

- AC adaptor
- Bluetooth® 2.0
- Bluetooth® BLE 4.0



350-10-8

Dual-ramp Platform

- 6.5 ft cable between scale base and indicator
- Dual ramp for easy access from either side
- Safety guide rails on two sides



350-10-7

Single-ramp Platform

- 6.5 ft cable between scale base and indicator
- Single ramp
- Safety guide rails on three sides

WHEELCHAIR SCALES

VACUUM FORMED SCALE TOP



Wheelchair Handrail Scales Built to enhance safety and support, Rice Lake's wheelchair scales with handrails simplify weighing for patients with limited mobility.



350-10-3 Dual-ramp Handrail

- Dual ramp
- Handrail for patient stability

All handrail models are designed to fold for easy portability and storage. Folded depth is 11.75in.



350-10-2 Single-ramp Handrail

- Single ramp
- Handrail for patient stability



350-10-4 Dual-ramp Handrail with Seat

- Includes padded, fold-away chair with back rest
- Dual ramp
- Handrail for patient stability



WHEELCHAIR SCALES

VACUUM FORMED SCALE TOP

Easy-to-Read Display

Each of Rice Lake's medical scales is equipped with a large, brightly lit indicator. The keys are strategically positioned on the overlay to aid operator ease while performing Body Mass Index (BMI) or tare entry. To transmit weight, height, and BMI data, the indicator comes equipped with a USB port or you have the option to upgrade to the optional Bluetooth and Bluetooth Low Energy (BLE) technologies.

The indicator can be powered by batteries or with an adapter, providing accurate weight information for both portable and stationary scales. The indicator is equipped with a battery life annunciator to help alert you when batteries need to be changed. When powered by batteries, the indicator will automatically dim the display back light, and you can adjust the auto-off feature to help extend battery life.

These indicators can also be configured to meet your specific needs, such as toggling the display between pounds and kilograms. You can also lock the display to measure in only pounds or kilograms at all times.



This 350-10-4 dual-ramp handrail with seat is folded for portability.

SALES AND TECHNICAL ASSISTANCE

RICE LAKE
WEIGHING SYSTEMS

HEADQUARTERS

230 West Coleman Street
Rice Lake, Wisconsin 54868 - USA
Tel: (715) 234 9171 | Fax: (715) 234 6967
www.ricelake.com